

A CUT ABOVE

INTRODUCING

Tri-Medics™

Precision Surgical Instruments

The moment you pick up a **Tri-Medics Precision Surgical Scissor** you feel the comfort, control, and smoothness that you are accustomed to with your precision instruments. You notice how the instrument is light, almost as if it is actually an extension of you, a part of your hand. When you cut with the **Tri-Medics Precision Surgical Scissor**, you experience the difference of an amazingly sharp, clean, accurate cut right to the tip.



*Ring Handle
Tenotomy Scissor*



*Spring Handle
Straight Scissor*



*Spring Handle
Curved Scissor*



*Spring Handle
Potts 45° Scissor*

- ▼ Consistently Sharp; Cuts at the Tip
- ▼ Lightweight
- ▼ Splay Resistant
- ▼ Durable, Strong, and Precise
- ▼ 100% Performance Guarantee

 **Tri-Medics™**
THE SHEAR ADVANTAGE


MADE IN THE
★ USA ★

Tri-Medics™ Patented Technology

Manufactured using a proprietary patented technology, **Tri-Medics Precision Surgical Instruments** are laser cut, stamped, and assembled from flexible, resilient stainless steel unlike conventional instruments which are forged from hardened, inflexible steel.

Flexible, resilient stainless steel provides two important advantages for surgeons. First, **Tri-Medics Precision Surgical Instruments** are lightweight and comfortable; many first-time users think the instruments are made of titanium. Second, **Tri-Medics Precision Surgical Scissors** consistently remain sharp and cut to the tip under normal use.

Comfort and performance; it's what makes a *precision surgical instrument* a **Tri-Medics**.



Straight Scissor



Curved Scissor



Potts 45° Scissor

Code	Description	Length
TMSHSS-7	Spring Handle Straight Scissor	7"
TMSHCS-7	Spring Handle Curved Scissor	7"
TMSHP45S-7	Spring Handle Potts 45° Scissor	7"
TMSHFNH-7C	Spring Handle Flat Needle Holder with Carbide	7"
TMSHFNH-7D	Spring Handle Flat Needle Holder with Diamond Dust	7"
TMSHFNH-8C	Spring Handle Flat Needle Holder with Carbide	8"
TMSHRNH-8C	Spring Handle Round Needle Holder with Carbide	8"
TMSHFNH-8D	Spring Handle Flat Needle Holder with Diamond Dust	8"
TMSHRNH-8D	Spring Handle Round Needle Holder with Diamond Dust	8"
TMRHS-95	Ring Handle Iris Scissor	95 mm
TMRHS-120	Ring Handle Tenotomy Scissor	120 mm
TMRHS-145	Ring Handle Stevens Tenotomy Scissor	145 mm

For more information on how you can have greater comfort and performance for your procedures, call your **Tri-Medics Representative at 888-490-5321**. To place an order, call **888-325-5551** or fax to **508-384-8424**.

**MADE IN THE
★ USA ★**



Tri-Medics and The Shear Advantage are trademarks of Tri-Medics LLC.
© 2011, Tri-Medics LLC. All Rights Reserved. P/N 001 Rev.02

Contact Info:

2 Hampshire Street
Suite 102
Foxboro, MA 02035
888-490-5321

www.tri-medics.com
info@tri-medics.com